

Lampiran 3. Cara pengujian konsentrasi TSS secara gravitimetri

1. Alat dan Bahan:

A. Alat:

- Kertas Saring
- Corong Kaca (8)
- Erlenmeyer 250 ml (8)
- Crush Tank (1)
- Pinset (1)
- Kaca arloji (1)
- Neraca analitik (1)
- Desikator (1)
- Oven (1)
- Stirer (1)
- Magnetic stirer (1)
- Gelas beaker 250 ml (4)

B. Bahan

- Aquades
- Sampel air limbah



